



ALMICRO[®] X10 : Multiviewing/Decahead Microscope

This Microscope comes with Infinity corrected plan optical system, which provides excellent optical quality image makes it easier to simultaneously see the slide and study it together for 10 people. This makes the microscope optimum for Medical/Dental College and Research use.

Specifications

- ✳ 1 Trinocular main Head & 9 Seidentopf Binocular heads inclined at 30°
- ✳ Wide Field Eyepiece EWF 10X/22mm
- ✳ Quadruple Nosepiece
- ✳ Infinity plan Objective 4X, 10X, 40X, 100X
- ✳ Abbe Condenser N.A.1.25
- ✳ Coaxial Coarse and Fine Adjustment
- ✳ Double Layer Mechanical Stage 180mmx160mm
- ✳ Halogen lamp 12V/100W, Brightness Adjustable
- ✳ LED Cursor Pointer, Brightness Adjustable